
LePurseal[®]

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1.General	
1.1 Company Profile	37
1.2 User Technical Safety Information	38
1.3 Liability disclaim	39
2.Introduction	
2.1 Safety Sign	40
2.2 Safety Instructions	41
2.3 Notice	42
2.4 Error description	43
2.5 Account permissions	44
2.6 Product introduction	45
2.6.1Parameter	45
2.6.2 Performance	46
3 Start-up	
3.1 Accessories list	49
4.Operation	
4.1 Operating procedures	51
4.2 Menu	52
4.3 Introduction to main interface icons	53
4.4 Emergency stop	54
4.5 Properly insert the tubing	55
4.6 Touch screen	56

4.6.1Self-test	56
4.6.2Login	58
4.6.3Sealing	58
4.6.4Export Data	61
4.6.5User management	63
4.6.6Settings	65
4.6.7Change password	66
4.6.8 Temperature Validation	67
4.6.9Device Information	68
4.7Equipment Cleaning.....	68
4.7.1Cleaning instructions	68
4.7.2Recommended cleaning	69
4.7.3Process description	69
5.Maintenance	
5.1Maintenance	70

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 A yellow triangular warning sign with a black border. Inside the triangle, there are three wavy lines representing heat or steam, and a horizontal line below them.	<p>当心高温表面，小心烫伤。封管的加热片表面在工作过程中存在高温，不得随意触摸。</p>
 A yellow triangular warning sign with a black border. Inside the triangle, there is a black silhouette of a hand being crushed between two vertical lines, representing a hand injury hazard.	<p>当心夹手，注意手部安全，避免发生手部意外伤害。</p>
 A yellow triangular warning sign with a black border. Inside the triangle, there is a large black exclamation mark, representing a general warning or hazard.	<p>本说明以及装置带有“警示标志”的“警告”是指，如果不遵守相关要求，不采取相应措施，就存在有可能造成死亡或严重人身伤害的潜在危险。</p>

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界面图片	描述
 The image shows the main splash screen of the LEPURE machine. It features the brand name 'LEPURE' in large white letters, with '全自动无菌封管机 II' (Fully Automatic Sterile Capping Machine II) below it. The background is dark blue with a glowing grid pattern and some light blue particles.	step1 插入电源插头，打开电源开关。
 The image shows a screen with a dark blue background and a glowing blue border. At the top, it says 'LEPURE' and '08:26:19'. The main text reads: '开始自检前 请移除夹管区域内所有物品 并合上保护盖!' (Before starting self-check, please remove all items from the clamping area and close the protective cover!).	step2 防护盖自动弹出，出现左侧界面。 *本机通电开机后，如果长时间不进行自检，本机将会进入保护状态，需断电重启封管机。
 The image shows the same screen as step 2, but with a blue button labeled '开始自检' (Start Self-check) at the bottom center.	step3 合上防护盖，界面出现“开始自检”按钮。

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1.General

1.1 Company Profile

Shanghai LePure Biotech Co., Ltd., established in 2011, is a high-tech enterprise engaged in the research and development, production, and sales of disposable consumables and equipment for the biopharmaceutical industry.

As one of the earliest domestic enterprises providing biopharmaceutical process solutions, LePure Biotech pioneered the localization of disposable solutions in the Chinese biopharmaceutical industry. For more than a decade, LePure Biotech has been committed to innovation in the industrial supply chain and high-quality domestic alternatives, always putting customers first and striving for excellence. LePure Biotech continuously develops new disposable bioprocess application products with a focus on quality. LePure Biotech holds multiple patents and technologies related to disposable bioprocess products for biopharmaceutical processes, covering upstream cell culture, downstream filtration, formulation filling, and fluid transfer. These products have multiple advantages, including safety, reliability, easy verification, low cost, and environmental friendliness, which help ensure the stability of drug processes and drug quality for biotechnology and pharmaceutical companies.

LePure Biotech is driven by technological innovation, fully leveraging the advantages of disposable technology in terms of efficiency, safety, convenience, and economy,

empowering biopharmaceuticals, providing high-quality process solutions, and committed to becoming the most trusted supplier for global biopharmaceutical companies.

1.2 User Technical Safety Information

This manual contains all the necessary information for the designated use of the product described in this document and is intended for technically qualified personnel. Qualified personnel are those authorized by the responsible person of system safety to perform necessary activities, who have the appropriate education, experience, and training, understand relevant standards, regulations, accident prevention ordinances, and operating conditions, and can identify and avoid hazards that may arise during these activities. The manufacturer is not liable for any damage or subsequent damage to the equipment or applications caused by accidents.

1.3 Liability disclaimer

The information provided in this manual is subject to change without prior notice. All statements and information in this manual are accurate and complete, but we do not provide any warranties or implied warranties for them.

Regardless of any other warranties made in this manual, all documents, files, and software in this manual may still have defects. We do not make any warranties, expressed or implied, including (but not limited to) warranties of merchantability, fitness for a particular purpose, non-infringement, and warranties arising from trade, use, or customary practices.

In no event shall we be liable for any indirect, special, consequential, or incidental damages, including (but not limited to) any loss of profits, data, or other losses resulting from the use or inability to use this manual. The user assumes full responsibility for the use of the product.

This equipment is intended solely for the purposes specified in this manual, and any unauthorized use or modification of this equipment without our prior written consent is considered improper operation. We shall not be liable for any damages or losses resulting from such actions, and the risk is assumed by the operator. We are not responsible for any damages or subsequent losses to the equipment or applications caused by accidents.

2.Introduction

2.1 Safety Sign

 A yellow triangular warning sign with a black border, containing three wavy lines representing heat.	<p>Be careful of the hot surface and avoid scalding. The heating surface of the sealing tube has high temperature during operation and should not be touched without authorization.</p>
 A yellow triangular warning sign with a black border, containing a hand being pinched by a mechanical part.	<p>Be careful not to pinch your hands and pay attention to hand safety to avoid accidental hand injuries.</p>
 A yellow triangular warning sign with a black border, containing a large exclamation mark.	<p>This notice and the "warning" with a "warning sign" on the device indicate that if you do not comply with the relevant requirements and take corresponding measures, there is a potential danger of death or serious personal injury.</p>

2.2 Safety Instructions

For the personal safety of operators and the safety use of the equipment, please read the following information carefully.

- Before turning on the machine, check if the rated voltage of the equipment is consistent with the local voltage.
- This device is not suitable for explosive environments.
- Not following the manufacturer's instructions may result in unreliability, dangerous situations, and safety interferences.
- This device is developed for laboratory use and must prevent any type of fluid or steam.
- The device can operate within the range of +18°C to +40°C. The sealing parameters are only applicable to room temperature conditions.
- The optimal relative humidity for sealing is $\leq 60\%$.
- Ensure that the equipment is not in a closed state during the sealing process.
- To improve process reliability, ensure that there is an uninterrupted power supply (UPS) available.
- Ensure that no foreign objects enter the internal components of the equipment during use to prevent damage.
- Make sure to log out after use, otherwise another user may use your account to operate the equipment.
- Not following the instructions in the manual may result in unreliable, hazardous situations, and safety interference.

2.3 Notice

- The power supply used for the equipment must be grounded.
- The equipment should be kept clean and dry. If it is dirty, it should be cleaned promptly to avoid affecting the normal operation of the equipment.
- When the equipment malfunctions, professional maintenance personnel should be notified to inspect and repair it.
- Interruptions during the sealing process can lead to issues and must be avoided. It is essential to strictly follow the operating instructions.
- The equipment must be disinfected before being sent to the supplier or manufacturer.

2.4 Error description

Code	Description	Solution
ERR1	Electric motor motion mechanism malfunction.	If any abnormal situation occurs during the operation of the device, please immediately turn off the power and restart it. If the problem persists, please contact our engineers for support and assistance.
ERR2	Anomaly in the fixed-end heating mechanism.	
ERR3	Abnormality in the mobile-end heating mechanism.	
ERR4	The fixed-end heating mechanism is not heating up.	
ERR5	The mobile-end heating mechanism is not heating up.	
ERR6	Cooling timeout alarm (cooling system abnormal).	

2.5 Account permissions

	ADMIN	Operator	Viewer	Example
Sealing	<input checked="" type="checkbox"/>	<input checked="" type="checkbox"/>	<input checked="" type="checkbox"/>	<input checked="" type="checkbox"/>
Exporting data	<input checked="" type="checkbox"/>	<input type="checkbox"/>	<input checked="" type="checkbox"/>	<input type="checkbox"/>
User management	<input checked="" type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>
Settings	<input checked="" type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>
Change Password	<input checked="" type="checkbox"/>	<input checked="" type="checkbox"/>	<input checked="" type="checkbox"/>	<input type="checkbox"/>
Device Information	<input checked="" type="checkbox"/>	<input checked="" type="checkbox"/>	<input checked="" type="checkbox"/>	<input checked="" type="checkbox"/>
Temperature Validation	<input checked="" type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>

* "Example" is a one-click login feature that allows you to log in without a password.

2.6 Product introduction

2.6.1 Parameter

Power	100~240VAC, 50/60Hz
System power	120W
Heat sealing time	Sealing tube for 3~5.5 minutes
Applicable pipeline specifications and models	<p>Thermoplastic tubing/BPT with an outer diameter of 6.4 mm ~ 25.4 mm (1/4" ~ 1") that can be sealed with or without liquid</p> <p>*The steam sterilized tubing should be left in place for 24 hours before use in case of any undesirable appearance caused by sealing.</p>
Dimension	Length * Width * Height: 37.5*16*15.5 cm
Net weight	5.5 kg
Operating environment	Temperature: +18°C ~ 40°C; Humidity ≤ 60%RH
Material	Base: 304 stainless steel, Housing: ABS engineering plastic, Flip cover: PMMA

Noise	≤65dB
Housing treatment	High-temperature baking paint 、 resistant to disinfectant wiping
Permission management	Three-level authority management

2.6.2 Performance

(1) Working principle

LePure Aseptic Sealing Machine has two corresponding heating plates. The tubing is placed between the two heating plates, clamped and heated. Under high temperature and pressure, the inner wall of the tubing is fused. After cooling, the tubing is completely sealed in the middle section. It can be cut with scissors at the cut line position in the middle when needed.

(2) Applicability

LePure Biotech aseptic tube sealing machine is suitable for heat sealing thermoplastic tubing with an outer diameter of 6.4 mm-25.4 mm. The heat sealing process does not require additional consumables and disposable technical products, such as disposable bioreactors and disposable storage bags, to be used in conjunction, which is economical and efficient.

(3) Temperature setting and tubing selection procedure

LePure aseptic sealing machine can seal different brands and materials of thermoplastic tubing, and adapt to different tubing materials through different temperature settings. The machine provides programs for multiple brands and different sizes of tubing. When ordering, we will set the corresponding parameters according to the confirmed tubing material to adapt to it.

The correspondence between brands and tubing models is detailed in the table

below:

Le-Flex®	
1/8"-1/4"	3/16"-5/16"
1/4"-3/8"	1/4"-7/16"
1/4"-1/2"	3/8"-5/8"
1/2"-3/4"	3/4"-1"

Advanta Flex®	
1/8"-1/4"	1/4"-3/8"
1/4"-7/16"	3/8"-5/8"
1/2"-3/4"	5/8"-7/8"
3/4"-1"	

C-Flex®	
1/8"-1/4"	1/4"-3/8"
1/4"-7/16"	3/8"-5/8"
1/2"-3/4"	3/4"-1"

BPT
BPT

Defoamer
Defoamer

*The steam sterilized tubing should be left in place for 24 hours before use in case of any undesirable appearance caused by sealing.

3 Start-up

3.1 Accessories list

Accessory name	Quantity	Option
Automatic Aseptic Sealer II	1	<input checked="" type="checkbox"/>

Power cord	1	<input checked="" type="checkbox"/>
Product manual	1	<input checked="" type="checkbox"/>
Power adapter	1	<input checked="" type="checkbox"/>

4.Operation

4.1 Operating procedures


-  This device has a heat source, please pay attention to the temperature of each part to avoid burns!
-  Keep the clamping chamber clean and free of contaminants, clean it regularly.If there is contamination on the heating plate, it will affect the

4.2 Menu


4.3 Introduction to main interface icons

	Tube sealing
	Export data
	User management
	Settings
	Change password
	Device Information

 A circular icon with a blue glow, containing a stylized graph with a temperature symbol (°C) and a checkmark.	Temperature Validation
---	------------------------

4.4 Emergency stop



4.5 Properly insert the tubing


*Please insert the tubing according to the following orientations:

4.6 Touch screen


4.6.1 Self-test

Interface picture	Describe
 The image shows the main logo screen of the LEPURE Automatic Aseptic Sealer II. The background is dark blue with a grid of glowing blue lines and dots, suggesting a technical or digital theme. The text 'LEPURE' is prominently displayed in a large, white, sans-serif font, with 'Automatic Aseptic Sealer II' written below it in a smaller, white font.	<p>step1</p> <p>Plug in the power supply and turn on the power switch.</p>
 The image shows a touch screen interface with a dark blue background. At the top, 'LEPURE' is on the left and '08:26:19' is on the right. A central white box with a blue border contains the text: 'Before starting the build-in self-test Please remove all items in the clamping area and close the protective cover!'. The text is centered and uses a white, sans-serif font.	<p>step2</p> <p>The protective cover pops up automatically, and the left interface appears.</p>
 The image shows the same touch screen interface as in step 2, but with an additional blue button at the bottom center labeled 'Start Self-test' in white text. The rest of the screen content remains the same.	<p>*After the machine is powered on, if it does not perform a self-test for a long time, it will enter a protection state and require powering off and restarting the tube sealing machine.</p> <p>step3</p> <p>Close the protective cover, and the 'Start Self-Check' button will appear on the interface.</p>


step4

After pressing the 'Start Self-check' button, the device begins to self-check.

step5

Self-check completed, left interface prompts.


step6

The login interface pops up, allowing access to the corresponding account for operation.

4.6.2 Login

Interface picture

Describe


After entering the login password, you can log in to the corresponding level and perform operations.

*ADMIN: Initial password is ADMIN; Example: One-click login; Other permissions need to be created by ADMIN.

4.6.3 Sealing


Interface picture

Describe


step1

Click on the 'Seal Tube' entrance on the main interface to enter the left interface.


step2

Please enter the batch number and select the brand and tube type from the drop-down menu.

step3

After confirming the brand and tube type, click 'Next' to display the left interface.

step4


Place the tube in the middle position of the tube holder. After covering the protective cover, click 'Start clamping'.

*Please place the tube flat and adhere it to the black tube placement platform before sealing the tube.

step5


Pre-clamp before sealing the tube.

step6


Tips before sealing the tubes.

*While on this screen, Please confirm that both sides of the tube have been clamped.


step7

Start preheating.


step8

Start sealing the tubes.

step9

Start cooling.

step10


Sealing completed.

- ✧ To seal with a different brand or tube type, select 'Select Tube Type' again.
- ✧ To seal with the same brand and tube type, select 'Continue Sealing'.


4.6.4 Export Data

Interface picture


Describe


1) Click on the 'Export Data' entrance on the main interface to enter the left interface. You can view and export sealing records and operation logs.


2) Click on 'Export Sealing Records' to enter the left interface. Select 'Previous Page' and 'Next Page' to switch between the previous and next page lists. Click on the input box in the middle of the bottom to quickly search for the specified page list.


3) The export options include the following filter selections; click the filter dropdown box, Select the 'All' option to export all the data.


4) Click on the filtering dropdown menu and select the "Batch Number" option to filter by batch number.


5) Click on the filtering dropdown menu and select the "Personnel" option to filter by personnel.

6) Click on the filtering dropdown menu and select "Time" option to filter by time.

7) Click on "Export Operation Log" to access the left-hand interface. Use the "Previous Page" and "Next Page" buttons to switch between pages of the list. The input box at the bottom center can be used to quickly search for a specific page in the list.


8) Click on the "Export" button in the upper right corner of step 7 to access the left-hand interface. Insert the USB flash drive and click on the "Export Data" button to export the data.


4.6.5 User management

Interface picture


Describe

1) Click on the "User Management" entry on the main interface to access the left-hand interface. Users with ADMIN privileges can access this interface and have the authority to create or delete users.


2) ADMIN is unique and cannot be deleted or modified.

3) Click on the "+ Create User" button in the upper right corner of Figure "2)" to enter the left interface. Enter the corresponding information to create a new user.

*The account number is unique. Accounts that have been used cannot be recreated.


4) Click on the management list to enter the left interface. You can modify and delete users from there.


4.6.6 Settings

Interface picture


Describe

1) Click the "Settings" button on the main interface to enter the left-hand interface. You can set the time, language, and tube type.

2) Time setting.

3) Language setting.


4) Management type setting.

4.6.7 Change password


Interface picture

Describe


1) Clicking on the "Change Password" button on the main interface will take you to the left-hand side interface. You can select your account and then proceed to change your password.

Password rules: 1-18 characters of any type

2) Password change success prompt.

3) Password change error prompt.


4.6.8 Temperature Validation

Interface picture	Describe
	<p>step1</p> <p>Click the "Temperature Validation" entry on the main interface to enter the left interface.</p>
	<p>step2</p> <p>Close the protective cover, and the 'start Temp. Validation' button will appear on the interface.</p> <p>*Only administrators have the authority to operate, and it is recommended that it be completed by a Lepure validation engineer.</p>
	<p>step3</p> <p>Press the 'Start Temp. validation' button to perform temperature validation.</p>
	<p>step4</p> <p>After the temperature validation is completed, the left interface will be displayed.</p>

4.6.9 Device Information

Interface picture

Describe

Click on the 'Device Information' entry to enter the left interface, where you can view some device information.

4.7 Equipment cleaning

4.7.1 Cleaning instructions

- The Automatic Aseptic Sealer II may not be cleaned with a spray- or steam procedure. Spray- or steam cleaning can cause a reduced life time of single components or lead to damage of single components.
- Clean the equipment using a lint-free cloth dipped in alcohol for disinfection and cleaning (refer to Figure 1 and Figure 2)

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4.7.2 Recommended Items

- Ethanol 70% (e.g. Klericide 70/30)
- Sodium hypochlorite 6% (e.g. Klericide-CR)

4.7.3 Process description

step1

The equipment self-check is complete. Keep the black cover open to clean the entire sealed area.

step2

"Wipe the sealing position with a clean tissue soaked in a cleaning agent (preferably non-woven fabric is recommended).

5.Maintenance

5.1Maintenance

- The equipment should be cleaned weekly, and the outer surface should be wiped with a dust-free cloth sprayed with alcohol before each use.
- Regularly check the sealing area for dirt and clean it immediately if any. Clean the sealing area with a lint-free cloth dipped in an appropriate amount of alcohol.
- The equipment's operation status should be checked regularly, and if any abnormalities are found, please contact the local sales and after-sales
- Perform annual equipment verification and maintenance.
- If the equipment needs to be returned to the factory for repair, please use our dedicated packaging to avoid transportation damage. If the unpacking package is lost, please contact the local sales representative for purchase.
- After-sales service hotline: 400-902-0506.

*Please contact your sales engineer for device maintenance and upkeep plans.

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