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Safety Information

Read all Installation and Operation Instructions before use, use Midcap according to the following instructions.

Warning

- Reduce the risk of pollutant exposure
- Use personal protection during installing and operating in accordance with standard operating procedures Open the inlet & outlet vales and release all the system pressure


Important Note

- Please follow the operation guide; otherwise, the operation will be invalid
- Do not exceed the operating parameters of the product
- Do not use online steam sterilization
- Midcap is disposable, not recommended for reuse

Procedures

- 1.Connect the inlet & outlet with TC25 interface and gasket, select a proper tube and peristaltic pump, close the outlet
- 2.Open the air vent of Midcap (Luer cap)
- 3.Start the pump, pump the purified water into the Midcap and exhaust the air inside
- 4.Pause the pump, close the air vent
- 5.Open the outlet, start the pump to flush Midcap (Flush Midcap according to the Flushing steps bleed)
- 6.After finishing the Flushing, pause the pump, the purified water or buffer is replaced into the liquid to be filtered, start the pump. (The flow rate is usually 50~ 60 ml/min depending on the liquid viscosity, the type and proportion of impurity)

7. Some liquid will remain in the filter after the filtration, which can be recovered with proper buffer or low-pressure clean air. The recommended air operating pressure: 0.3~0.5 bar; Max. pressure < 1 bar; recovery time: 5~10 min

Flushing

- Usually 1.5 L of purified water is used for flushing Midcap before use. Recommend flushing flow rate: 75~150 mL/min; Max. inlet pressure < 1 bar; Flushing time: 10~20 min
- Flush according to the correct flow direction
- It is recommended to flush with buffer after flushing with purified water. The operation process and parameters are the same as those of purified water

Technical Parameters

Parameter list of LeClafine® Midcap

Effective Filtration Area	300 cm ²
Inlet&Outlet	TC25
Max. Inlet Pressure	3.5 bar(51 psi)
Max. Operate Temp	40°C (104° F)
Max. Steam sterilization Temp/Cycle	126°C (259° F)/ 1 cycle
Max. ΔP	2.4 bar(35 psi)
Housing	PP / Glass fiber

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