

**LEPURE**

# Extractables Guide

**LFRS Platinum Cured Silicone Tubing**

## Catalogue

1.Introductions .....	3
2.Purpose and Methodology of Extractables Testing.....	5
3.Information of Component and Instruments .....	6
3.1 Information of Component .....	6
3.2 Information of Instruments .....	6
4.USP < 665 > COMPLIANCE.....	7
4.1 Overview of Extractables Protocol .....	7
4.2 Results of Extractables Testing.....	8
4.2.1 Results of Elemental Impurities .....	9
4.2.2 Results of Organic Compounds .....	10
5.Safety Assessment .....	11
Reference .....	17

## 1. Introductions

The pharmaceutical single use system produced by Shanghai LePure is widely used in the bio-pharmaceutical process. Its main application fields include the research, development and production of antibody, vaccine and cell therapy products. At present, pharmaceutical single use systems are widely used in upstream, downstream and final filling processes. Therefore, end users must fully understand and verify their interactions with bio-pharmaceutical solutions and final drug products. In order to ensure the product quality, the pharmaceutical enterprises shall conduct comprehensive analysis and testing in the early stage of process development and within the scope of process monitoring and quality control, so as to prove the purity, efficacy and safety of drugs. Product safety assessment shall be carried out in the process of design and development of the single use system. These validation studies involve complete quantitative, identification and toxicological assessment of the leachables, which are the substance remained in the drug solution due to the interaction between the drug solution and the pharmaceutical single use system. Leachables are a subset of the extractables that can be extracted from pharmaceutical single use systems. Usually, the solvents and conditions for extractables testing are more severe than that for leachables. The purpose of this Guide is to provide the worst-case data on extractables to support the validation studies conducted by process developers and toxicologists.

The safety problem of disposable components, which is the most common problem in the pharmaceutical single use system, has always been the most concerned problem of many pharmaceutical enterprises, especially for the safety of extractables. For this kind of disposable components, LePure has referenced the technical data of foreign raw material suppliers, the technical guidelines for compatibility studies published by domestic CDE (including Technical Guidelines for Compatibility Studies of Chemical Injection and Pharmaceutical Glass Packaging Container (Tentative), Technical Guidelines for Compatibility Studies of Chemical Injection and Plastic Packaging Material (Tentative), and Technical Guidelines for Compatibility Studies of Chemical Injection and Elastomeric Seals (Tentative) in China; and the Application and Technical Guidelines for Single Use System (Tentative) published by the China Center for Food and Drug International Exchange in November 2017), the technical guidelines formulated by relevant SUS organizations, and USP <665> and USP <1665> and developed a reasonable test scheme for the dissolved substances in single use system.

The potential dissolved substances in the pharmaceutical single use system produced may come from surfactants, lubricants and additives in the process of plastic processing, or from the shedding of raw materials of material structure and oligomer monomer. This study report summarizes the information on extractables from disposable components in which includes elements and organic compounds (nonvolatile compounds, semi-volatile and volatile compounds, small-molecule volatiles) in three kinds of simulated solvents, mainly referring to USP<665>, BPOG guide for pharmaceutical single use system. Elements were detected by inductively coupled plasma-mass spectrometry (ICP-MS), non-volatile compounds were detected by Ultra-High performance liquid chromatography (HPLC), semi-volatile and volatile compounds were detected by gas

chromatography-mass spectrometry (GC-MS), and small-molecule volatiles were detected by headspace gas chromatography-mass spectrometry (HS-GC-MS).

LEPURE Biotech

## 2. Purpose and Methodology of Extractables Testing

According to the guidelines of USP < 665 > and USP < 1665 >, the extractables test scheme was formulated respectively, and the extractables were tested according to this scheme. Before the test, LePure has confirmed the analysis and evaluation threshold (AET). According to the guidelines of Product Quality Research Institute (PQRI), AET was determined as a threshold. When the concentration of a compound exceeds the threshold, the compound should be identified, quantified and reported. In addition, the toxicity for this compound needs to be evaluated. AET is obtained through conversion according to the appropriate safety assessment threshold (SCT) or toxicological concern threshold (TTC), taking into consideration the dosage of the product. When the concentration of a compound is lower than the threshold, it can be considered that the toxicity of the compound is very low and will not be harmful to human. As LePure single use systems may be exposed to a variety of drugs and chemical reagents, and the maximum daily dosage of drugs cannot be confirmed at this stage, based on the minimum detection limit of the instrument, and according to the guide of BPOG, the report limit of 0.1µg/mL was defined as AET for organic compounds, and the report limit of 20ng/mL was set as AET for inorganic substances for this study. Converted to the concentration for surface area, AET for organic compounds was 0.016µg/cm<sup>2</sup>, and AET for inorganic substances was 0.003µg/cm<sup>2</sup>.

### 3. Information of Component and Instruments

#### 3.1 Information of Component

Table 1 Information of Component

Name	Material	Cat. No.	Batch. No.
LFRS Platinum cured silicone tubing	Silicon	LFRS2416	20201203

Note: This guide is applicable to other components constructed from the same materials.

#### 3.2 Information of Instruments

Table 2 Information of Instruments


Name	Model/Specification	Manufacturer
Gas Chromatography-Mass Spectrometry	8890-5977B	Agilent
Ultra-High Performance Liquid Chromatography	Vanquish Flex	Thermo Fisher
Inductively Coupled Plasma-Mass Spectrometry	iCAP RQ	Thermo Fisher

## 4. USP < 665 > Compliance

### 4.1 Overview of Extractables Protocol

USP<665> guide was referred to determine the extraction solution, extraction method and test instrument used in extractable study. An overview of extractables study process is provided in Table 3.

Table 3 Extractable Study Process



## 4.2 Results of Extractables Testing

The elemental, non-volatile, semi-volatile and volatile extracts in LFRS platinum cured silicone tubing were detected. In addition, we focused on the antioxidants and their degradation products, fatty acids, phthalate plasticizers, polycyclic aromatic hydrocarbons, lubricants, siloxanes, vulcanizing agents, nitrosamines and other additives during the experiment and data analysis.

LEPURE Biotech



## 4.2.1 Results of Elemental Impurities

Table 4 Results of Elemental Impurities

Element	ICH Q3d Class	Conc. ( $\mu\text{g}/\text{cm}^2$ )	
		pH3	pH10
Cd	Class 1	<LOR	<LOR
Pb	Class 1	<LOR	<LOR
As	Class 1	<LOR	<LOR
Hg	Class 1	<LOR	<LOR
Co	Class 2A	<LOR	<LOR
V	Class 2A	<LOR	<LOR
Ni	Class 2A	<LOR	<LOR
Pt	Class 2B	<LOR	<LOR
Li	Class 3	<LOR	<LOR
Sb	Class 3	<LOR	<LOR
Ba	Class 3	<LOR	<LOR
Mo	Class 3	<LOR	<LOR
Cu	Class 3	<LOR	<LOR
Sn	Class 3	<LOR	<LOR
Cr	Class 3	<LOR	<LOR
B	N/A	0.006	0.005
Sr	N/A	<LOR	<LOR
W	N/A	<LOR	<LOR
Mg	N/A	0.004	<LOR
Al	N/A	<LOR	0.004
Ca	N/A	0.012	0.008
Ti	N/A	<LOR	<LOR
Mn	N/A	<LOR	<LOR
Fe	N/A	0.004	<LOR
Zn	N/A	<LOR	<LOR
Si	N/A	1.583	1.915

Note: "LOR" means limit of report, and the concentration of LOR is  $0.003\mu\text{g}/\text{cm}^2$ .

## 4.2.2 Results of Organic Compounds

Table 5 Results of Organic Compounds

Model Solvent	Analytical Method	Compound Name	CAS	Conc. ( $\mu\text{g}/\text{cm}^2$ )
50% EtOH	GC-MS	Unknown	N/A	0.089
	HS-GC-MS	Octamethyl cyclotetrasiloxane	556-67-2	0.194
pH3	N/A	ND	ND	ND
pH10	N/A	ND	ND	ND

Note: 1) "LOR" means limit of report and the concentration of LOR is  $0.016\mu\text{g}/\text{cm}^2$ .

2) "ND" means not detected.

## 5. Safety Assessment

The toxicity of extractables and leachables must be evaluated for the effects on both patients and process. Although almost any quantity of certain compounds in a drug is considered unacceptable (e.g., ICH Q3C class-1 solvents), the toxicity of extractables or leachables must be observed in the broader context of the following criteria, the actual concentration of the leachables in the final drug, the mode of administration, the dosage, the duration of treatment, the number of patients, and the risk benefit assessment.

Therefore, the toxicity is not only related to the identification and concentration of the extractables, or only related to the amount of leachables in the process fluid or pharmaceutical intermediates. The daily intake of patients can be obtained taking into consideration information of extractable concentration, model solvent volume, test component contact area, process component contact area, process batch and dosage. The daily intake of single compound can be compared with PDE value. For compounds or unknown substances whose PDE value cannot be obtained, the worst scenario can be assumed.

- 1) All extractables are migrating to the final product.
- 2) All extractables are considered as DNA reactive impurities (genotoxicant).

The purpose of determining the toxicological concern threshold (TTC) in ICH M7 is to define a common acceptable exposure level for compounds that have gone through toxicological studies (see Table 6). An appropriate limit value can be chosen based on treatment cycle and treatment route to conduct safety assessment.

Table 6 Acceptable Intakes for an Individual Impurity

Duration of treatment	≤1 month	>1-12 months	>1-10 years	>10 years to lifetime
Daily intake [µg/day]	120	20	10	1.5

Depending on the toxicity categorization and concentration of each detected compound and elements, there is no high-risk compounds or elements detected for LFRS platinum cured silicone tubing. The toxicological assessment of extractables or leachables for drug products should be performed based on the process conditions and clinical dose for patient.

## Appendix 1: Study of Elements Impurity

Table 7 Instrument Method for Elements Impurity by ICP-MS

Mass Spectrometry Conditions	Parameter	Numerical Value					
	Plasma Flow	14 L/min					
	RF Power	1550 W					
	Detection Model	STD/KED					
	Solution Rinse Time	60s					
	Solution Lift Time	60s					
	Target Element	Element	Mass	Element	Mass	Element	Mass
		Li	7	Cu	63	Ni	60
		B	11	Zn	66	As	75
		Mg	24	Sn	118	Sr	88
	Al	27	Ba	137	Mo	95	
	Ca	44	Si	29	Cd	111	
	Ti	48	V	51	Sb	121	
	Cr	52	Mn	55	W	182	
	Fe	57	Co	59	Pt	195	
	Hg	202	Pb	208			
Internal Standard Element	Element	Mass	Element	Mass	Element	Mass	
	Be	9	Y	89	Lu	175	
	Sc	45	Rh	103	Bi	209	
	Ge	73	In	115			

## Appendix 2: Study of Organic Compounds

Table 8 Scan Method for Semi-Volatile Compounds by GC-MS

Column	DB-5HT (30 m *0.25mm, 0.10µm)		
Injection Volume	1 µL		
Carrier Gas	He		
Temperature	Injection Temp.: 250°C, Transferline Temp.: 280°C, MS Source Temp.: 230°C		
Flow Control Mode	Constant flow, rate is 1.0 mL/min		
Acquisition Type	Full Scan 33-650		
<b>Temperature Program</b>			
No.	Rate (°C/min)	Temp. (°C)	Hold Time (min)
1	/	50	1
2	10	160	6
3	5	315	11

Table 9 Scan Method for Volatile Compounds by HS-GC-MS

Headspace Condition	80°C 30min		
Column	DB-624UI, 30m×0.25mm, 1.4µm		
Injection Volume	1mL		
Carrier Gas	He		
Split Flow Rate	1mL/min		
Temperature	Injection Temp.: 250°C, Transferline Temp.: 250°C, MS Source Temp.: 280°C		
Flow Control Mode	Split ratio: 20:1		
Acquisition Type	MS, Full Scan: 30-500		
<b>Temperature Program</b>			
No.	Rate (°C/min)	Temp. (°C)	Hold Time (min)
1	/	50	1
2	7	250	30

## Appendix 3: Study of Specially Concerned Compounds

Table 10 Instrument Method for Phthalate Plasticizers by GC-MS

Column	DB-5HT (30 m *0.25mm, 0.10 $\mu$ m)		
Injection Volume	1 $\mu$ L		
Carrier Gas	He		
Temperature	Injection Temp.: 250°C, Transferline Temp.: 280°C, MS Source Temp.: 230°C		
Flow Control Mode	Constant flow, rate is 1.0 mL/min		
Acquisition Type	<b>SIM</b>		
	<b>Target Compounds</b>	<b>m/z</b>	
	Naphthalene	127, 128, 129	
	Acenaphthylene, DMP, Acenaphthene, Fluorene, DEP	149, 152, 153, 163, 165, 178	
	Phenanthrene, Anthracene, DIBP	149, 176, 178, 223	
	DBP, DMEP, Fluoranthene, Pyrene, BMPP, DEPP, DPP,	45, 59, 149, 193, 202	
	BBP, DHXP, Benzo[a]anthracene, <b>Chrysene</b> , DBEP, DCHP, Diphenyl phthalate, DEHP	104, 149, 225, 228	
	Benzo[b]fluoranthene, Benzo[k]fluoranthene, DNOP, Benzo[a]pyrene, DNP	57, 149, 252	
	Indenopyrene, Dibenzoanthracene, Benzoperylene	138, 276, 278	
	<b>Temperature Program</b>		
<b>No.</b>	<b>Rate (°C/min)</b>	<b>Temp. (°C)</b>	<b>Hold Time (min)</b>
1	/	50	1
2	10	160	6
3	5	315	11

Table 11 Instrument Method for Other Non-volatile Compounds by UPLC-DAD

Column	ACQUITY UPLC BEH C18 (100mm×2.1mm, 1.7µm)		
Injection Volume	5µL		
Mobile Phase A	Ultra-pure water		
Mobile Phase B	Acetonitrile		
Flow rate	0.3mL/min		
Temperature	Sample tray: 5°C, column: 40°C		
UV Wavelength	230nm and full wavelength scanning		
<b>Gradient</b>			
Gradient	Time (min)	Mobile Phase A (%)	Mobile Phase B (%)
	0.0	95	5
	6.0	95	5
	20.0	0	100
	30.0	0	100
	30.1	95	5
	40.0	95	5



## Reference

1. USP<665> Plastic Components and Systems Used to Manufacture Pharmaceutical Drug Products and Biopharmaceutical Drug Substances and Products
2. USP<1665> Characterization and Qualification of Plastic Components and Systems Used to Manufacture Pharmaceutical Drug Products and Biopharmaceutical Drug Substances and Products
3. BioPhorum Best Practices Guide for Extractables Testing of Polymeric Single-Use Components Used in Biopharmaceutical Manufacturing
4. BioPhorum Best Practices Guide for Evaluating Leachables Risk from Polymeric Single-Use Systems Used in Biopharmaceutical Manufacturing
5. ICH Q3D(R1) Guidelines for Elemental Impurities
6. ICH M7(R1) Evaluate and Control DNA Reactive (Mutagenic) Impurities in Drugs to Limit Potential Carcinogenic Risk
7. Application and Technical Guide of Disposable Use System.

# LEPURE

LePure Biotech

## LEPURE

Shanghai LePure Biotech Co., Ltd.

Website: [www.lepure-bio.com](http://www.lepure-bio.com)

Building 3, 410 Yunzhen Road, Songjiang, Shanghai, China 201600

E-mail: [marketing@lepure-bio.com](mailto:marketing@lepure-bio.com)