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## Introduction

The following procedures should be followed for the installation and use of LePure LeSiever® Sterilizing Capsule Filter. These instructions and the information contained within the product document includes brochure, datasheet, validation guide must be read thoroughly. It is critical to follow these instructions to ensure optimal filter performance. If some of the procedures do not suit your needs, please consult LePure Biotech. LePure Biotech shall not be held accountable for any damage to the capsule filters caused by improper operation by any person.

## Specification

Please check the brochure, chemical compatibility chart, datasheet and validation guide for details, or contact LePure Biotech.

## Sterilization

Caution: Steam-in-place [SIP] is not applicable for capsule filters

- The autoclave cycle should ensure that the filters are sterilized and must be validated accordingly.
- The valve protective caps, and Tri-Clamp protective cap of the capsule filter should be removed before autoclaving.
- The vent and drain valves should be opened before autoclaving.
- Sterilization pouches ensure seal integrity and circulation of hot steam.
- Capsules should be oriented in a direction that flow indicator arrows are pointing down during autoclave cycle.

**The LeSiever® Sterilizing capsule Filter should only be sterilized as follows:**

Autoclave	125°C ,30 mins,5 cycles
-----------	-------------------------

## Notes:

1. After autoclave sterilization, the filter can only be used until the temperature falls below 60° C.
2. Sterilization using process parameters outside the validation range should be independently validated.

## Installation

LeSiever® capsules must be oriented in the process so that process flow follows the flow indicated by the flow direction arrow on the capsule, and the pressure gauge is connected to the production pipeline according to the connection mode of the inlet and outlet of the capsule filter and the nozzle size for real-time detection of pressures.

## Filter Wetting

Filters should be installed in an appropriate orientation, which follows the flow indicated by the flow direction arrow, close the outlet valve of the capsule filter, slowly inject liquid (WFI, buffer or material solution) into the filter chamber at a pressure of 2 psi. Open the upper vent valve until the liquid flows out of the vent valve, then close the vent valve. Open the downstream drain valve of capsule, and adjust water flow to approximately 2 LPM/10 inch of filtration area for a minimum 5 minutes.

## Filter Flushing

When the filter is subjected to autoclaving, flushing after autoclaving is highly recommended. Flushing is beneficial to filtrate quality and will reduce leachables from the filter. Recommended Flush Volume (Liters) per LeSiever® Filters:

Model of Capsule Filter	Minimum flushing volume
Size 1	0.4 L
Size 5	1 L
2 Inch	2 L
4 Inch	4 L
5 Inch	5 L
10 Inch	10 L
20 Inch	20 L

## Integrity Testing

Sterilizing Capsule Filter should be integrity tested post-use, according to a GMP requirement for many applications. Pre-use, post-sterilization integrity testing is recommended to ensure that the filter is capable of performing its stated function. Contact LePure Biotech for recommended integrity test procedures and integrity test values.

## Capsule Filter Replacement

Filter capsules should be replaced in line with the GMP requirements of the process. LeSiever® capsules are disposable. No attempt should be made to clean disposable filter capsules. Where filter capsules are used for more than one manufacturing batch, replacements are recommended when the maximum allowable differential Capsule Filter pressure has been reached, if the flow rate has become unacceptable or if the cumulative steam life has been reached, whichever occurs first. Discard Capsule Filter in accordance with local environmental procedures. The maximum allowable differential pressure is shown in the table below:

	Size1	Size5/2"/4"/5"/10"/20"
Maximum upward pressure differential	4.0 bar	5.0 bar
Maximum downward pressure differential	2 bar	2.5 bar

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United States LePure Biotech LLC